

CARPOFIX MANU

Fixed wrist hand brace

Ref. SRT 203

Indications

It is used in sprains, dislocations, arthrosis, carpal tunnel syndrome, radial nerve paralysis, tenosynovitis, fractures at the level of the area, both as orthopedic and postoperative treatment.

Description

Brace made of textile, provided with:

- one mouldable metal palm splint
- 2 inserted plastic splints on the dorsal side of the hand.


Hook and loop fastener with 3 straps provided with buckles, on the dorsal side of the hand.

Length: 24 cm.


The orthosis stabilizes and immobilizes the radiocarpal joint.

Composition

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



air transfer



water transfer



Product details



- hook and loop fasteners




- mouldable metal splint



- Available color:
black

Available sizes, left/right versions

Sizes	1	2	3	4	5
Wrist circumference (cm)	15-17	17-19	19-21	21-23	23-25
Reference	SRT203M1	SRT203M2	SRT203M3	SRT203M4	SRT203M5



After the Reference must be added **L** for Left, **R** for Right

The circumference of the wrist is measured