


HALLUFIX

Hallux valgus corrector

Ref. SRT 301

INDICATIONS

It is used for the correction of hallux valgus and post surgery.
To be worn during night.



MEDICAL




DAILY

DESCRIPTION


Orthosis with anatomical design, made of textile, provided with a mouldable metal splint designed to physiologically correct the position of the hallux. Hook and loop fastener.

COMPOSITION

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



AIR FLOW



WATER
TRANSFER



AVAILABLE COLOURS

Black



AVAILABLE SIZES

Universal size, for left and right versions.
Forefoot with circumference between 20 – 29 cm.

