

TRIAGEN

Mobile knee brace with lateral articulations

Ref. SRT 319

Indications

It is used in lateral instability of the knee, traumas or injuries at the level of the knee, inflammatory processes, for postoperative recovery.

Description

The orthosis is made of textile, open in the front, provided with:

- a non-slip silicone elastic band (10mm) inside the upper part
- 2 metal splints with mobile articulations, one on each side
- 2 straps provided with elastic hook and loop fastener
- cutout in the popliteal region for comfort.

Hook and loop fastening on the front side.


Length: 32 cm.

The orthosis is a support for the knee joint, the 2 polycentric, lateral metal splints providing stability and safety in use.

Due to its particular shape, it can to be used also in case of a deformed knee. The compression can be adjusted by the frontal closure system and the 2 fastening straps.

Composition

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



air flow

water transfer



- Available colour: black

Available sizes

Sizes	1	2	3	4	5	6	7	8
Knee circumference (cm)	33-35	35-37	37-39	39-41	41-43	43-45	45-47	47-49
Reference	SRT319M1	SRT319M2	SRT319M3	SRT319M4	SRT319M5	SRT319M6	SRT319M7	SRT319M8

The circumference of the knee is measured at the middle of the knee cap



Product details



- orthosis open in the front
- non-slip silicone elastic band



- the popliteal area is free



- metal splints with mobile joint

