


PATELLAFIX 01

Patellar tendon support

Ref. SRT 307

INDICATIONS


It is used to reduce discomfort and inflammation by compression at the level of the patellar tendon.



MEDICAL



SPORT



DAILY

DESCRIPTION


Orthosis with anatomical design, made of textile, provided with:

- one pad that exert compression at the level of the patellar tendon
 - posterior protective sleeve, optional.
- Hook and loop back fastener.


Width: 6 cm.

COMPOSITION

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



AIR FLOW



WATER
TRANSFER

AVAILABLE COLOURS

Black



AVAILABLE SIZES

Universal size

For knee circumference between 33 – 49 cm.

