

AKTIV LIGTH 01 TRIA

Mobile knee brace with patellar support

Ref. SRT 308

Indications

It is used in cricks, luxations, post consolidated fractures, inflammatory processes, lateral knee joint instability. It is recommended both as orthopedic treatment and postoperative.

Description

It is a tubular, closed orthosis, made of textile, provided with:

- a non-slip silicone elastic band (10mm) inside the upper part
- 4 spring type metal splints, 2 on each side
- 2 straps for fixation, provided with elastic hook and loop fastener. The elastic fastening follows the conical shape of the thigh, preventing the orthosis from slipping
- patellar support (polyethylene foam).

Length: 29 cm.

The orthosis is a support for the knee joint, the 2 spring type metal splints for each side providing stability and safety in use.

Due to its tubular, sleeve shape, it is not recommended to be used in case of oedema or significant swelling.

Composition

Neopren coated on both sides with fabric, 100% polyester.




Product details



- tubular, closed orthosis
- hook and loop elastic fastener



- patellar support made of polyethylene foam



- spring type metal splints



- Available colour: black

Available sizes

Sizes	1	2	3	4	5	6	7	8
Knee circumference (cm)	33-35	35-37	37-39	39-41	41-43	43-45	45-47	47-49
Reference	SRT308M1	SRT308M2	SRT308M3	SRT308M4	SRT308M5	SRT308M6	SRT308M7	SRT308M8

The circumference of the knee is measured at the middle of the knee cap

