


GENU AKTIV

Mobile knee orthosis

Ref. SRT 349

INDICATIONS

It is used in case of light injuries of the knee, inflammatory processes, sports injuries.



MEDICAL



SPORT



DAILY

DESCRIPTION

Orthosis with anatomical design, open on the back and made of textile, provided with:

- patellar support made by overlapping 2 layers of material
- 3 fixing straps.


Hook and loop fastener.

Length: 25 cm.


The orthosis stabilizes the knee joint and supports the patella, without exerting pressure on it.

COMPOSITION

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



AIR FLOW



WATER
TRANSFER



AVAILABLE COLOURS

Black



AVAILABLE SIZES

Universal size

For knee circumference between 33 – 45 cm.

