


ARMFIX EASY

Upper limb support

Ref. SRT 239

INDICATION

It is used to support the arm and the forearm.



MEDICAL



SPORT



DAILY

DESCRIPTION

Support with anatomical design, made of one single textile piece that supports the arm, provided with:


- one adjustable support strap (with "crocodile" type fastener)
- one fix strap to support the wrist.

Hook and loop fastener.


The orthosis is an effective support for the arm and the forearm.

COMPOSITION

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



AIR FLOW



WATER
TRANSFER

AVAILABLE COLOURS

Black



AVAILABLE SIZEZ

Universal size

For body height between 150-190 cm.



REMODELED

