


ELBISPORT

Mobile elbow brace

Ref. SRT 241

INDICATIONS

It is used in sprains, luxations and contusions.



MEDICAL



SPORT



DAILY

DESCRIPTION


Orthosis with anatomical design, made of textile.
Hook and loop fastener.

Length: 18 cm.


The orthosis is a support for the elbow, stabilizes and protects it.

COMPOSITION

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward



AIR FLOW



WATER
TRANSFER



AVAILABLE COLOURS

Black



AVAILABLE SIZES

Universal size

For elbow circumference between 22 – 34 cm.

