

## ELBI PREMIUM

Epicondylitis brace

Ref. SRT 207

### INSTRUCTIONS FOR USE



Dear Customer!

Thank you for the confidence you have shown in purchasing Triamed products! We would like to inform you that Triamed products are made of the best performing materials, all having quality certificates and biocompatibility tests. These orthoses comply with the Technical Specification and the requirements of Regulation (EU) MDR 2017/745 of the European Parliament on medical devices and comply with the requirements set out in Annex I - General Requirements for safety and performance- and Annex VIII - Classification rules. Conformity is declared on the basis of Article 52 (7), Annex II and Annex III of MDR 2017/745.

#### Important information

The orthosis is used according to the recommendations of the specialized physician.

Choosing the size - For optimal results, the appropriate size will be selected based on the measurements table of the product. A size too small may lead to intolerance to wearing or the damage of the product! A size too large significantly diminishes the efficacy of the orthosis.

This orthosis should only be applied in accordance with the instructions mentioned here.

The product is intended for single patient use.

No liability is accepted for incorrect product use.

Guarantee is bound to statutory regulations.

#### Available sizes

Universal size

For elbow circumference between 22 – 34 cm.

#### Indications

It is used in epicondylitis, tendinitis, elbow joint oversteering (tennis elbow).

#### Description

Orthosis with anatomical design, made of textile, provided with:

- pressing pad at the epicondyle, that exerts pressure on the ligaments
- fastening strap with buckle.

Hook and loop fastener.

Width: 8 cm.

Due to the fastening strap with buckle the orthosis can be adjusted to the optimal compression degree, so that the pressing pad exerts the necessary pressure on the ligaments.

#### Composition

Neoprene coated on both sides with fabric:

PES Polyester - 100%; SBR Styren-butadien Rubber – 100%; PA Nylon – 82% + EA Spandex – 18%.

#### Positioning

The orthosis is placed on the upper third of the arm, the band is passed through the buckle and it is adjusted to the optimal compression.

#### Cleaning



To be washed by hand with warm water, no bleaching, no spinning, no ironing, no chemical cleaning!

#### Storage

At room temperature, in a dry space with good ventilation.

#### Commercial warranty of the product: 30 days.

**Date of production:** see packaging label.

**Validity date:** no expiry date.

To validate the warranty, the seller must indicate the date of purchase. The buyer is obliged to keep and present the warranty coupon.

In order to solve the problems that appeared in the warranty, the following data must be notified: Article reference, batch and date of manufacture, indicated on the packaging box of the product.

#### Deadline for communication of the solution: 15 days.

Please send notifications and complaints to the following email address: [suportclienti@triamed.ro](mailto:suportclienti@triamed.ro)

<b>WARRANTY COUPON</b>	
30 days from the date of purchase	Batch and date of manufacture
Product reference	Date of purchase
The reason of complaint	Seller name and signature



**Manufacturer**

**SC TRIAMED SRL**

**Production** Calea Buziaşului nr. 77, Hall 1-C1/A | Timișoara | Timiș | România

**Registered office** Uliuc 262C | Sacoșu Turcesc | Timiș | România

**Phone** 00 40 256 48 73 88

office@triamed.ro | triamed@gmail.com | www.triamed.ro

Edition 4, Revision 4/31.01.2025