

GENU MOBIL

Mobile knee brace with patellar support

Ref. SRT 305

Indications

It is used in case of light knee injuries, inflammatory processes.

Description

It is a tubular orthosis made of textile, provided with patellar support (polyethylene foam).

Length: 29 cm.

The orthosis stabilizes the knee joint and supports the knee cap without compressing it.

Composition

Neoprene coated on both sides with fabric, 100% polyester.



Product details



- tubular, closed orthosis
- patellar support made of polyethylene foam



- Available color: black

Available sizes

Sizes	1	2	3	4	5	6	7	8
Knee circumference (cm)	33-35	35-37	37-39	39-41	41-43	43-45	45-47	47-49
Reference	SRT305M1	SRT305M2	SRT305M3	SRT305M4	SRT305M5	SRT305M6	SRT305M7	SRT305M8

The circumference of the knee is measured at the middle of the knee cap

