


## ELBISPORT AKTIV

Mobile elbow brace with elastic strap

Ref. SRT 245

### INDICATIONS

It is used in sprains, luxations and contusions.



MEDICAL



SPORT



DAILY

### DESCRIPTION


Orthosis with anatomical design, made of textile, provided with an elastic band for adjusting the compression. Hook and loop fastener.

Length: 18 cm.


The orthosis is a support for the elbow, stabilizes and protects it.

### COMPOSITION

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



AIR FLOW



WATER  
TRANSFER



### AVAILABLE COLOURS

Black



### AVAILABLE SIZES

Universal size

For elbow circumference 22 - 34 cm.



## Available sizes

Sizes	1	2	3	4	5	6
Arm circumference (cm)	22-26	26-30	30-34	34-38	38-42	42-46
Reference	<b>SRT238M1</b>	<b>SRT238M2</b>	<b>SRT238M3</b>	<b>SRT238M4</b>	<b>SRT238M5</b>	<b>SRT238M6</b>



The circumference of the arm is measured at the proximal 1/3