

SINTEX JUNIOR

Arm support

Ref. SRT 224

Indications

It is used to support the arm and the forearm.

Description


The support is made of one single textile piece that supports the arm, provided with 2 adjustable support straps.

Hook and loop fastener.

The orthosis is an efficient support for the arm and forearm. Due to the two fastening straps, it does not allow the arm abduction and external rotation.

Composition

The rigid textile material, made of 100% polyester, is a spatialized 3D mesh, similar to a honeycomb, allowing a very high air permeability in both directions.



air transfer



also available in adult version



Product details



- provided with adjustable support straps
- support strap for the thumb




- spatialized 3D mesh



- Available color: black

Available sizes

Sizes	1	2	3
Forearm length (cm)	16-18	18-20	20-22
Reference	SRT224M1	SRT224M2	SRT224M3



The length of the forearm is measured