

## AKTIV

Mobile knee brace with patellar support and lateral articulations

Ref. SRT 304

### Indications

It is used in lateral instability of the knee, traumas or injuries at the level of the knee, inflammatory processes, for postoperative recovery.

### Description

It is a tubular, closed orthosis, made of textile, provided with:

- a non-slip silicone elastic band (10mm) inside the upper part
- 2 metal splints, one on each side
- 2 straps of elastic fastening material. The elastic fastening follows the conical shape of the thigh, preventing the orthosis from slipping
- patellar support (polyethylene foam).

Length: 32 cm.

The orthosis is a support for the knee joint, the 2 polycentric, lateral metal splints providing stability and safety in use.

Due to its tubular, sleeve shape, it is not recommended to be used in case of oedema or significant swelling.

### Composition

Neoprene coated on both sides with fabric, 100% polyester.



### Product details



- tubular, closed orthosis
- hook and loop elastic fastener



- patellar support made of polyethylene foam



- metal splints with mobile joint



- Available colour: black

### Available sizes

Sizes	1	2	3	4	5	6	7	8
Knee circumference (cm)	33-35	35-37	37-39	39-41	41-43	43-45	45-47	47-49
Reference	SRT304M1	SRT304M2	SRT304M3	SRT304M4	SRT304M5	SRT304M6	SRT304M7	SRT304M8

The circumference of the knee is measured at the middle of the knee cap

