

## GENU MOBIL

Mobile knee brace with patellar support

Ref. SRT 305

### Indications

It is used in case of ligh knee injuries, inflammatory processes.

### Description

It is a tubular orthosis made of textile, provided with a non-slip silicone elastic band (10mm) inside the upper part and patellar support (polyethylene foam).

Length: 29 cm.

The orthosis stabilizes the knee joint and supports the knee cap without compressing it.

### Composition

Neoprene coated on both sides with fabric, 100% polyester.



### Product details



- tubular, closed orthosis
- patellar support made of polyethylene foam



- Available colour: black

### Available sizes

| Sizes                   | 1        | 2        | 3        | 4        | 5        | 6        | 7        | 8        |
|-------------------------|----------|----------|----------|----------|----------|----------|----------|----------|
| Knee circumference (cm) | 33-35    | 35-37    | 37-39    | 39-41    | 41-43    | 43-45    | 45-47    | 47-49    |
| Reference               | SRT305M1 | SRT305M2 | SRT305M3 | SRT305M4 | SRT305M5 | SRT305M6 | SRT305M7 | SRT305M8 |

The circumference of the knee is measured at the middle of the knee cap

