


## TRIAFIX ARM 03

Arm support

Ref. SRT 247

### INDICATIONS

It is used to support the arm and the forearm.



MEDICAL



SPORT



DAILY


### DESCRIPTION

Support with anatomical design, made of one single textile piece that supports the arm, provided with adjustable support strap (with "crocodile" type fastener and protection sleeve). The length of the textile piece which supports the arm is adjustable by folding the support at the wrist, on the marked stitches. After measuring the length of the forearm (olecranon-radiocarp), the support is folded as needed. Hook and loop fastener.


The orthosis is an effective support for the arm and the forearm.

### COMPOSITION

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



AIR FLOW



WATER  
TRANSFER



### AVAILABLE COLOURS

Black



### AVAILABLE SIZES

Universal size

Forearm length between 22-40 cm.

