


## HIPSPORT

Hip orthosis

Ref. SRT 353

### INDICATIONS

It is used to protect the hip joint, in case of minor muscle injuries, during sustained sports activities, minor traumas.



MEDICAL



SPORT



DAILY

### DESCRIPTION

Orthosis with anatomical design, open in the front, made of textile material:


- one adjustable component for the hips (waist)
- one adjustable component for the thigh.

Hook and loop fastening on the front at the level hip and at the thigh level, with 2 straps.


The waist component is positioned with the help of 2 straps with hook and loop over the thigh component, depending on hip circumference  
Length: 39 cm.

### COMPOSITION

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



AIR FLOW



WATER  
TRANSFER

### AVAILABLE COLOURS

Black



### AVAILABLE SIZES

Universal size

For hip circumference between 75 – 125 cm.

