


## ABD SPORT

Thoraco abdominal support

Ref. SRT 150

### INDICATIONS

It is used in case of abdominal injuries, postoperative or after childbirth.



MEDICAL



SPORT



DAILY


### DESCRIPTION

Support with anatomical design, made of elastic fabric. Hook and loop fastener on the front side. Width: 24 cm.


It is a traditional elastic corset that, due to the material, allows normal thoraco abdominal breathing, while ensuring uniform compression on the entire covered surface.

### COMPOSITION

Elastic fabric: 40%polyester, 40% polyamide, 20% elastomer, with a structure that allows high air permeability in both directions.



AIR FLOW



WATER  
TRANSFER



### AVAILABLE COLOURS

Black



### AVAILABLE SIZES

Universal size

For waist circumference between 65 – 125 cm.

