


ANKLEFIX ACTIV

Mobile ankle support with plastic side splints

Ref. SRT 314

INDICATIONS

It is used in cricks, luxations, post consolidated fractures, inflammatory processes, after the removal of the orthopedic plaster cast.



MEDICAL



SPORT



DAILY

DESCRIPTION

Orthosis with anatomical design, made of 2 rigid, plastic side splints, with protection pads made of polyethylene foam, covered with soft textile and is provided with:

- 2 lateral fastening straps provided with hook and loop
- one fastening strap at heel level that allows the adjustment of the width.


Height: 25 cm.

The orthosis is a support for the ankle, the 2 rigid, plastic side splints ensure stability and safety in use.


COMPOSITION

Side splints made of polypropylene.

The protection pads are manufactured by thermoforming and have a shape that follows the anatomical line of the ankle.



AIR FLOW



WATER
TRANSFER

AVAILABLE COLOURS

Black



AVAILABLE SIZES

Universal size, for left and right versions

For ankle with circumference between 20 – 30 cm.

