

THORATEX 03

Elastic thoracic support

Ref. SRT 114

Indications

It is used to correct vicious positions of the thoracic spine or in case of kyphosis.


Description

The support is made of elastic fabric, provided with:


- 2 paravertebral integrated plastic rods
 - 2 elastic straps which surround the shoulders and are fixed on the chest, correcting the position of the trunk by traction. The straps are provided with sleeves for underarm protection.
- Hook and loop fastener on the front side.

Composition

Elastic fabric: 50% polyester, 33% polyamide, 17% elastane, with a structure that allows high air permeability in both directions.



air transfer




also available in Junior version



Product details



- hook and loop fastening on the front side
- sleeves for underarm protection



- elastic fabric with a structure that allows high air permeability in both directions



- Available colors: black, beige

Available sizes

Sizes	1	2	3	4	5	6
Torso circumference (cm)	75-85	85-95	95-105	105-115	115-125	125-135
Reference	SRT114M1	SRT114M2	SRT114M3	SRT114M4	SRT114M5	SRT114M6

The circumference of the torso is measured

