

ABD JUNIOR 01

Support for abduction

Code: **SRT 322**

Indications

It is used in the case of congenital hip dysplasia.

Description

It is made of textile, provided with a rigid pad for keeping legs in the abduction.

Hook and loop front closure.

Composition

Textile fabric - 100% cotton.


Rigid pad - 100% polyethylene foam.

Properties of the fabric


The cotton fabric is soft and does not cause skin irritation.

It allows good breathing to the skin and keeps the body dry.

It easily absorbs and removes body moisture.



airflow



water transfer

Available sizes

Sizes	1	2	3	4	5
Weight of the child (kg)	2-3	3-5	5-7	7-9	9-11
Code	SRT322M1	SRT322M2	SRT322M3	SRT322M4	SRT322M5

The weight of the child is measured