


## GENU MOBIL

Mobile knee orthosis with lateral metal articulations

Ref. SRT 350

### INDICATIONS

It is used in lateral knee joint instability, injuries, cricks, luxations, inflammatory processes. It is recommended both as orthopedic treatment and postoperative.



MEDICAL



SPORT



DAILY

### DESCRIPTION

Orthosis with anatomical design, open at the back, made of textile material, provided with:


- 2 metal splints with mobile articulation, one on each side
- patellar support made by overlapping 2 layers of material.

Hook and loop fastening on the front and on the back side. Length: 30 cm.


The orthosis is a support for the knee joint, the 2 metal splints providing stability and safety in use. It can be used also in case of a deformed knee, due to the possibility of adjusting the compression through the rear closing mode.

### COMPOSITION

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



AIR FLOW



WATER  
TRANSFER

### AVAILABLE COLOURS

Black



### AVAILABLE SIZES

Universal size

For knee circumference between 33 – 45 cm.

