

TRIAFIX ELBOW

Fixed elbow wrist hand brace

Ref. SRT 214

Indications

It is used in tendinopathy, carpal tunnel syndrome, radial nerve paralysis, cubital tunnel syndrome, fracture at the level of the wrist or elbow, both as orthopedic treatment and postoperative.

Description

Brace made of textile, provided with:

- one mouldable metal palm splint
- 2 fixed metal elbow splints set at 90°, mouldable, of the same length as the orthosis
- 3 straps with hook and loop (2 on the forearm and one on the arm)
- adjustable support strap (with "crocodile" type fastener and protection sleeve), optional. The positioning is done on the first fixing strap on the forearm and over the opposite shoulder.


Hook and loop fastener on the dorsal side of the hand.

Length: 51 - 57 cm.


The orthosis stabilizes and immobilizes the elbow and the radiocarpal joint.

Composition

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



air transfer



water transfer



- Available color: black

Available sizes, left/right versions

Sizes	1	2	3	4	5	6
Elbow circumference (cm)	22-24	24-26	26-28	28-30	30-32	32-34
Forearm length (cm)	22-25	22-25	25-28	25-28	28-31	28-31
Reference	SRT214M1	SRT214M2	SRT214M3	SRT214M4	SRT214M5	SRT214M6

After the Reference must be added **L for Left, R for Right**

The circumference of the elbow and the length of the forearm are measured



Product details



- fastening strap that limits the extension of the elbow



- fastening strap around the radiocarpal joint

