


UNOFIX SPORT

Fixed thumb brace

Ref. SRT 249

INDICATIONS

It is used in case of sprains, luxations and fractures of the thumb.



MEDICAL



SPORT




DAILY

DESCRIPTION


Brace with anatomical design, made of textile, the thumb is fixed with a plastic splint.
Hook and loop fastener.
Width: 10 cm.
The orthosis stabilizes and immobilizes the metacarpophalangeal and interphalangeal joint of the thumb.

COMPOSITION

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



AIR FLOW



WATER
TRANSFER

AVAILABLE COLOURS

Black



AVAILABLE SIZES

Universal size
For wrist circumference between 15 – 25 cm.

