


OMBIFIX SPORT

Abdominal support for umbilical hernia

Ref. SRT 149

INDICATIONS

It is used in case of abdominal disorders and umbilical hernia.



MEDICAL



SPORT



DAILY

DESCRIPTION

Support with anatomical design, made of elastic fabric, provided with one semirigid, standard dimension (polyethylene foam, Ø 11 cm), removable pad, on the front, at the level of the navel.

Hook and loop straps.

Width: 14 cm.


It is a traditional elastic corset that, due to the material, allows normal thoraco abdominal breathing, while ensuring uniform compression to on the entire covered surface.

The semirigid pad is detachable, it can be repositioned or removed. It provides protection and compression to the umbilical region. The support can be used in case of umbilical hernia, both before and after surgery.


COMPOSITION

Elastic fabric: 40% polyester, 40% polyamide, 20% elastomer, with a structure that allows high air permeability in both directions.

Abdominal pad: polyethylene foam lined with soft textile.



AIR FLOW



WATER
TRANSFER

AVAILABLE COLOURS

Black



AVAILABLE SIZES

Universal size

For waist circumference between 65 – 125 cm.

