

GENU REHA

Mobile knee brace

Ref. SRT 344

Indications

It is used to stabilize the knee, in case of ligament instability, genu recurvatum.

Description

The orthosis is made of a metal frame, provided with:

- 2 hook and loop rigid straps on the front
- rear band that allows the adjustment of the hyperextension blocking angle.

Composition

The metal frame is made of 3 mm duralumin.

The 3 rigid straps for fixation are made of rigid fabric, provided with protection pads.



Product details



- 2 rigid straps on the front




- rear band for the adjustment of the hyperextension blocking angle



- Available color: black

Available sizes

Sizes	1	2	3
Knee circumference (cm)	25-33	34-43	44-51
Reference	SRT344M1	SRT344M2	SRT344M3



The circumference of the knee is measured at the middle of the knee cap