


## CALFSPORT

Leg support

Ref. SRT 351

### INDICATIONS

It is used in leg muscle lesions, sport injuries.



MEDICAL



SPORT




DAILY

### DESCRIPTION


Orthosis with anatomical design, open on the front and made of textile.  
Hook and loop fastener on the front side with 3 straps.  
Length: 26 cm.

### COMPOZIȚIE

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moist to outward.



AIR FLOW



WATER  
TRANSFER



### AVAILABLE COLOURS

Black



### AVAILABLE SIZES

Universal size

For calf circumference between 29 – 48 cm.

