

TRIACARP M FINGER

Fixed wrist hand finger brace

Ref. SRT 213

Indications

It is used in sprains, dislocations, arthrosis, carpal tunnel syndrome, radial nerve paralysis, tenosynovitis, fractures at the level of the area, both as orthopedic and postoperative treatment.

Description

Brace made of textile, provided with:


- one mouldable metal palm splint
 - one mouldable metal thumb splint
 - strap around the radiocarpal joint, provided with hook and loop fastener.
- Hook and loop fasteners on the dorsal side of the hand.

Length: 24 cm.


The orthosis is a semirigid medical device that successfully stabilizes and immobilizes the radiocarpal joint, limiting flexion and extension, while the thumb metal splint immobilizes the thumb.

Composition

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



air transfer



water transfer



Product details



- strap around the wrist for a better fit
- hook and loop fastener




- mouldable metal thumb splint
- mouldable metal palm splint



- Available color: black

Available sizes, left/right versions

Sizes	1	2	3	4	5
Wrist circumference (cm)	15-17	17-19	19-21	21-23	23-25
Reference	SRT213M1	SRT213M2	SRT213M3	SRT213M4	SRT213M5



After the Reference must be added **L for Left, R for Right**

The circumference of the wrist is measured