

## ANKLEFIX CONFORT

Mobile ankle support with back fastening

Ref. SRT 311

### Indications

It is used in cricks, luxations, post fractures, after consolidated fractures, inflammatory processes, after the removal of the orthopedic plaster cast.

### Description

The orthosis is made of textile, provided with 2 fastening straps.

The heel area is open.


Hook and loop fastener on the back side.

Height: 20 cm.


The orthosis is a support for the ankle, ensuring stability while wearing.

### Composition

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



air transfer



water transfer



### Product details



- hook and loop fastener on the back side



- open heel area



- Available color: black

### Available sizes

Sizes	1	2	3	4	5
Ankle circumference (cm)	20-22	22-24	24-26	26-28	28-30
Reference	SRT311M1	SRT311M2	SRT311M3	SRT311M4	SRT311M5

The circumference of the ankle is measured around the malleolus

