

AKTIV LIGTH 02 TRIA

Mobile knee brace with patellar support

Ref. SRT 316

Indications

It is used in cricks, luxations, post consolidated fractures, inflammatory processes, lateral knee joint instability. It is recommended both as orthopedic treatment and postoperative.

Description

It is a tubular, closed orthosis, made of textile, provided with:

- a non-slip silicone elastic band (10mm) inside the upper part
- 2 spring type metal splints, one on each side
- patellar support (polyethylene foam).

Length: 29 cm.

The orthosis is a support for the knee joint, the two spring type metal splints for each side providing stability and safety in use.

Due to its tubular, sleeve shape, it is not recommended to be used in case of oedema or significant swelling.

Composition

Neoprene coated on both sides with fabric, 100% polyester.



Product details



- tubular, closed orthosis



- spring type metal splint



- Available colour:
black

Available sizes

Sizes	1	2	3	4	5	6	7	8
Knee circumference (cm)	33-35	35-37	37-39	39-41	41-43	43-45	45-47	47-49
Reference	SRT316M1	SRT316M2	SRT316M3	SRT316M4	SRT316M5	SRT316M6	SRT316M7	SRT316M8

The circumference of the knee is measured at the middle of the knee cap

