


## ARMSPORT

Shoulder orthosis

Ref. SRT 242

### INDICATIONS

It is used in light contusion, sprain, luxation, periarthritis humeroscapularis. The orthosis covers the shoulder and the arm and it is an elastic support for the shoulder joint and it can be used both as orthopedic treatment in minor disorders and for recovery in post immobilization period.



MEDICAL



SPORT



DAILY


### DESCRIPTION

Orthosis with anatomical design, made of textile material, provided with a fastening band under the opposite armpit, which has a protection sleeve. Due to the fastening strap provided with 2 buckles the orthosis can be both worn on the left and the right shoulder. The length can be adjusted by cutting at the level of the stitches to avoid fabric breaking, and repositioning the "crocodile" type fastener. Hook and loop fastener.


It is adjusted to the optimal compression, and can be used in sports activities that overload the shoulder joint.

### COMPOSITION

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



AIR FLOW



WATER  
TRANSFER

### AVAILABLE COLOURS

Black



### AVAILABLE SIZES

Universal size

For arm circumference between 22 - 46 cm.

**ONE  
SIZE**

