

## GEROMED 01

Mobile knee brace

Ref. SRT 318

### Indications

It is used in cricks, luxations, post consolidated fractures, lateral knee joint instability, inflammatory processes. It is recommended both as orthopedic treatment and postoperative.

### Description

The orthosis is open at the back, made of textile, provided with:

- a non-slip silicone elastic band (10mm) inside the upper part
- 4 spring type metal splints, 2 on each side
- patellar support (polyethylene foam).

Hook and loop fastening on the front side.


Length: 26 cm.

The orthosis is a support for the knee joint, the 4 spring type metal splints, 2 for each side providing stability and safety in use.


Due to its particular shape, it can be used also in case of a deformed knee. The compression can be adjusted by the back closure system.

### Composition

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



air flow



water transfer



- Available colour: black

### Available sizes

Sizes	1	2	3	4	5
Knee circumference (cm)	33-36	36-39	39-42	42-45	45-49
Reference	SRT318M1	SRT318M2	SRT318M3	SRT318M4	SRT318M5

The circumference of the knee is measured at the middle of the knee cap




### Product details



- orthosis open on the back side
- hook and loop fastener



- patellar support made of polyethylene foam



- spring type metal splints

